





## APPLICATION NOTE

High-resolution, live-imaging of *C. elegans* germline using *vivoChip*<sup>TM</sup>**Purpose**


The study of chromatin organization and defects in embryonic development of *C. elegans* using the *vivoChip*, which allows for high-resolution imaging of the entire 3-dimensional (3D) structures of the germline in *C. elegans*.

**Methods**

- Up to 40 adult animals are immobilized and oriented in the *vivoChip*-2x within 3 min.
- High-resolution images of the germlines are captured in 3D.
- The entire germline is visualized using two-color fluorescence reporters (see Fig 1).
- Chromatin organization/defects are identified and dead/viable embryo populations are quantified from the live time-lapse imaging of immobilized animals.



**Figure 1.** Images of AUM1039 animals using 20 $\times$ , 0.75NA objective in the Lionheart BioTek. The animals express a membrane bound GFP (*pie-1::gfp::PH*) and a nuclear mCherry (*pie-1::mCherry::his-58*)<sup>1</sup>.



**Figure 2.** AUM1039 germlines from two animals immobilized inside the parallel channels showing oocytes, spermatheca, and pachytene structures in two z-stack images.

**Figure 3.** Adult *C. elegans* in presence and absence of FUdR, an inhibitor for DNA synthesis. Bottom panel: Time-lapse images of the ROI (red box) collected at 2.1 s intervals show movement of viable embryos.

**Conclusions**

- The *vivoChip* enables uniform orientation of animals during immobilization for high-resolution imaging of entire germline organization.
- Anesthetic-free, time-lapse imaging provides a unique method to identify viable embryos.
- The *vivoChip* facilitates statistical evaluation of worm size, total number of embryos, and embryo viability for the development and reproductive toxicity (DART) assessment.
- Live imaging using *vivoChip* also allows researchers to investigate mechanisms of toxicity by identifying defects across different stages of meiosis and embryonic development.

<sup>1</sup>AUM1039 ItIs37 [P-*pie-1::mCherry::his-58* (pAA64) + *unc-19(+)*]; ItIs38 [P-*pie-1::GFP::PH* (PLC1delta1) + *unc-19(+)*] was gifted by Prof. Swathi Arur, MD Anderson.